



**TETRA**  
MEDICAL SUPPLY CORP.

6364 West Gross Point Road  
Niles, IL 60714  
800-621-4041  
www.tetramed.com

## Tetra-Flex Clip Free Woven Elastic Bandage

**Our Best Just Got Better!**



**Our Tetra-Flex bandage has been preferred by top orthopedic surgeons, head athletic trainers and key clinicians for years when compression and durability are an important component to treatment and therapy.**

**We have taken our finest quality, heavy-duty woven bandage and made it even better!**

**No more metal clips** - Tetra-FlexCF is a premium weave of spandex and selected long staple cotton providing maximum compression and durability, and it now has a **NEW SELF-CLOSURE SYSTEM!** Easy to apply and adjust. Durable to withstand repeat washings and applications.

**What's New:** We have reengineered our closure that is unique to our premium construction for maximum and effective clip free closure (core end) that firmly secures the bandage in place. Up to this point, only lighter weight knit bandages were available in clip free. Now you can have a high quality bandage where maximum compression and durability are a critical component in post-op procedures or treatment.

**Time Savings:** Our new clip free system makes it quick and easy to apply and adjust the bandages for both the clinician and the patient and eliminates the need for tape or other materials to keep it in place.

**Sizing:** Non-Sterile and Sterile in 5.5 yards (stretched length) in four standard widths (2", 3", 4" and 6") and in 11 yards (stretched length) in 4" and 6" widths.

**\*SEE ATTACHED FOR ORDERING INFORMATION.**

Tetra Medical Supply Corp. is a WBENC certified Women Owned Small Business and we welcome the opportunity to help organizations reach their diversity goals.



# TETRA

MEDICAL SUPPLY CORP.

6364 West Gross Point Road  
 Niles, IL 60714  
 800-621-4041  
 www.tetramed.com

Ordering Information for Tetra-Flex™CF (Clip Free) Woven Elastic Bandages.

**NOTE:** The new clip free bandages will start transitioning into inventory in early 2017.

Prev. Tetra #	New Tetra #	Description	UOM	Qty per UOM
0151-20	0151-2V	TETRA-FLEX CF,CLIP FREE WOVEN ELASTIC BANDAGE, LF*, 2" x 5.5 YD	DZ	12 rolls
0151-30	0151-3V	TETRA-FLEX CF,CLIP FREE WOVEN ELASTIC BANDAGE, LF, 3" x 5.5 YD	DZ	12 rolls
0151-40	0151-4V	TETRA-FLEX CF,CLIP FREE WOVEN ELASTIC BANDAGE, LF, 4" x 5.5 YD	DZ	12 rolls
0151-60	0151-6V	TETRA-FLEX CF,CLIP FREE WOVEN ELASTIC BANDAGE, LF, 6" x 5.5 YD	DZ	12 rolls
0155-24	0155-24V	TETRA-FLEX CF,CLIP FREE WOVEN ELASTIC BANDAGE, LF, 4" x 11 YD	DZ	12 rolls
0155-26	0155-26V	TETRA-FLEX CF,CLIP FREE WOVEN ELASTIC BANDAGE, LF, 6" x 11 YD	DZ	12 rolls
0151-2S	0151-2SV	<b>STERILE</b> TETRA-FLEX CF, CLIP FREE WOVEN ELASTIC BDG, LF, 2" x 5.5 YD	BX	10 rolls
0151-3S	0151-3SV	<b>STERILE</b> TETRA-FLEX CF, CLIP FREE WOVEN ELASTIC BDG, LF, 3" x 5.5 YD	BX	10 rolls
0151-4S	0151-4SV	<b>STERILE</b> TETRA-FLEX CF, CLIP FREE WOVEN ELASTIC BDG, LF, 4" x 5.5 YD	BX	10 rolls
0151-6S	0151-6SV	<b>STERILE</b> TETRA-FLEX CF, CLIP FREE WOVEN ELASTIC BDG, LF, 6" x 5.5 YD	BX	10 rolls
S155-24	S155-24V	<b>STERILE</b> TETRA-FLEX CF, CLIP FREE WOVEN ELASTIC BDG, LF, 4" x 11 YD	BX	10 rolls
S155-26	S155-26V	<b>STERILE</b> TETRA-FLEX CF, CLIP FREE WOVEN ELASTIC BDG, LF, 6" x 11 YD	BX	10 rolls

\*LF = Latex Free

Rev. 04/2017

Tetra Medical Supply Corp. is a WBENC certified Women Owned Small Business and we welcome the opportunity to help organizations reach their diversity goals.